



Single Patient Composite Disc Kit
3/8" Disc Kit Item Code: J75-1457
1/2" Disc Kit Item Code: J75-1458
Single Patient Composite Disc Kit helps reduce cross contamination and is conveniently packaged in individual, single patient units. Packaged 48/box, each unit contains four color coded discs and a disposable mandrel. Grits contained are one each of Coarse, Medium, Fine and Superfine. Price: **\$100.50**



Single-Use Diamond Composite Polishing Kit
Item Code: J75-1260
10 ea. cup, point, disc in 2 grits. Total - 60 pieces.
Price: **\$84.75**



Two-Step Diamond Composite Polishing Kit
Item Code: J75-1282
1 ea. sm point, point, cup, knife edge in 2 grits. Total - 8 pieces.
Price: **\$58.70**

Single-Use Diamond Composite Polishers

This disposable polishing system is ideal for the finishing and polishing of all composites, especially the industry's latest generation of composites. Polishers contain diamond abrasives in two different grits.

48 per Pkg.
Rec. RPM:
5,000-10,000



	Coarse			Fine		
Desc:	Disc	Cup	Point	Disc	Cup	Point
Item Code:	J75-1262	J75-1264	J75-1266	J75-1268	J75-1270	J75-1272
Shank:	RA	RA	RA	RA	RA	RA
Price:		\$66.75			\$99.75	

Two-Step Diamond Composite Polishers

Two-step diamond polishing system for finishing, smoothing and high-luster polishing of the latest generation of composites. Diamond abrasives are embedded in the silicone polishers.

3 per Pkg.
Rec. RPM:
5,000-10,000



	Coarse				Fine			
Desc:	Small Point	Point	Cup	Knife Edge	Small Point	Point	Cup	Knife Edge
Item Code:	J75-1284	J75-1286	J75-1288	J75-1290	J75-1292	J75-1294	J75-1296	J75-1298
Shank:	RA	RA	RA	RA	RA	RA	RA	RA
Price:		\$20.95				\$20.95		



Complete Disc/Strip Kit
Item Code: J75-1210

200 discs - 25 each 3/8" and 1/2" in Superfine, Fine, Medium and Coarse grits; 50 strips - 25 each Regular Coarse/Medium and Regular Fine/Superfine; 2 RA mandrels. Total - 252 pieces.
 Price: **\$79.50**

Composite Discs and Strips

Easy to use and ideal for all composites. Produce surface smoothness and visual gloss for more natural and life-like results. Color-coded to simplify your polishing procedures.

100 per Pkg.
 Rec. RPM:
 20,000

Desc:



Coarse



Medium



Fine



Superfine

Item Code:

J75-1212

J75-1214

Disc Refills 3/8"

J75-1216

J75-1218

Item Code:

J75-1220

J75-1222

Disc Refills 1/2"

J75-1224

J75-1226

Price:

\$38.95

Finishing and Polishing Strips

For interproximal trimming, finishing and polishing of composites. 100 per Pkg. 100 per pkg. Price **\$27.50**



Regular Coarse/Medium Strips / J75-1228

Narrow Coarse/Medium Strips / J75-1230




Regular Fine/Superfine Strips / J75-1232

Narrow Fine/Superfine Strips / J75-1234

Diamond Strips

6" long, extra fast cutting diamond strips available in superfine, fine & medium. 10 per Pkg. Price: **\$43.50**



Super Fine / J75-1236

Fine / J75-1240

Medium / J75-1244

Perforated design provides improved visibility, flexibility and control in reduction and contouring. 10 per Pkg. Price: **\$48.00**



Super Fine / J75-1238

Fine / J75-1242

Medium / J75-1245

Strip Saws

7" stainless steel saw for opening embrasures 8mm wide.



Desc: 8mm Strip Saw

Item Code: J75-1250

6/Pkg

Price:

\$20.25

Mandrels



Desc: RA Mandrel

Item Code: J75-1246

HP Mandrel

J75-1248

3/Pkg

3/Pkg

Price:

\$8.25

\$8.25

\$8.25

303 Mandrel

J75-1249

12/Pkg

\$10.50

Diamond Paste

Containing actual diamond particles for a superb, natural appearance. 2gm per pkg.



Desc: Extra Fine

Item Code: J75-0530

Fine

J75-0531

Regular

J75-0532

Price:

\$16.95

\$16.95

\$16.95

3 pack Asstd.

J75-0534

\$46.50

Separating Discs

7/8" X .009 - Made with the purest silicon carbide abrasive for fast cutting and finishing of porcelain interproximals and margins. Will not discolor porcelain. 25 per Pkg. Max. RPM: 30,000 **J75-1252** Price: **\$8.65**

7/8" X .015 - Made with the purest silicon carbide abrasive to provide fast cool cutting of porcelain and gold without discoloration. 25 per Pkg. Max. RPM: 30,000 **J75-1254** Price: **\$8.65**

7/8" X .025 - Ruby - Made with aluminum oxide that has a ceramic coating for long lasting, fast cutting of all chrome, precious and semi-precious metals. 100 per Pkg. Max. RPM: 30,000 **J75-1256** Price: **\$7.50**

7/8" X .025 - Stripper - Made with brown aluminum oxide for fast cutting of all chrome, precious and semi-precious metals. 100 per Pkg. Max. RPM: 30,000 **J75-1258** Price: **\$11.95**



Amalgam/Gold Polishing Kit
Item Code: J75-1300
 1 ea. mini point, point, cup, bullet in
 3 grits for polishing amalgam and
 metal. Total-12 pieces Price: **\$17.25**







Diamond Porcelain Polishing Kit
Item Code: J75-1326
 1 ea. sm point, point, cup, disc
 in 3 grits for polishing of
 ceramic, metal and zirconia.
 Total-12 pieces. Price: **\$112.50**



Standard Porcelain Polishing Kit
Item Code: J75-1358
 1 ea. sm point, point, cup, disc
 in 3 grits for polishing of
 porcelain and enamel.
 Total-12 pieces. Price: **\$19.50**






Amalgam/Gold Polishers

Three-step polishing system for removing, smoothing and polishing amalgam and gold. The silicon carbide is embedded in a medium-soft, flexible synthetic rubber matrix.

12 per Pkg. Rec. RPM: 7,000-12,000				
Desc:	Mini Point	Point	Cup	Bullet
Item Code:	J75-1302	J75-1304	J75-1306	J75-1308
Item Code:	J75-1310	J75-1312	J75-1314	J75-1316
Item Code:	J75-1318	J75-1320	J75-1322	J75-1324
Shank:	RA	RA	RA	RA
Price:			\$11.65	






Diamond Porcelain Polishers

Three-step diamond polishing system for removing, finishing and high-luster polishing of ceramic, metal and zirconia. Diamond abrasive is embedded in a hard-flexible, silicone polisher. The high-luster polish is achieved without paste.

3 per Pkg. Rec. RPM: 7,000-12,000					
Desc:	Wheel	Small Point	Point	Cup	Disc
Item Code:	J75-1328	J75-1330	J75-1332	J75-1334	J75-1336
Item Code:	J75-1338	J75-1340	J75-1342	J75-1344	J75-1346
Item Code:	J75-1348	J75-1350	J75-1352	J75-1354	J75-1356
Shank:	HP	RA	RA	RA	RA
Price:	\$39.00			\$22.50	

Standard Porcelain Polishers

Silicone polishers made specifically for polishing porcelain and enamel, with three grits for pre polishing, polishing and high shine super polishing.

10 per Pkg. Rec. RPM: 7,000-12,000					
Desc:	Wheel	Small Point	Point	Cup	Disc
Item Code:	J75-1360	J75-1362	J75-1364	J75-1366	J75-1368
Item Code:	J75-1370	J75-1372	J75-1374	J75-1376	J75-1378
Item Code:	J75-1380	J75-1382	J75-1384	J75-1386	J75-1388
Shank:	HP	RA	RA	RA	RA
Price:			\$11.25		

Ceramic Abrasives Kit
Item Code: J75-1448
 1 ea. disc, cone in
 2 grits.
 Total - 4 pieces.
 Price: **\$37.50**



Silicone Acrylic Polishing Kit
Item Code J75-1390
 1 ea. sm point, point, large point,
 jumbo in 3 grits. Total-12 pieces.
 Price: **\$44.95**



Ceramic Abrasives

High performance synthetically bound diamond abrasives for removing metal, specially designed for zirconia and alumina restorations without generating heat. Coarse abrasives for gross reduction, medium abrasives for finishing. For polishing, use Porcelain Polishers.

3 per Pkg.
 Rec. RPM:
 8,000-12,000



Coarse



Medium



Desc:
 Item Code:
 Shank:
 Price/3 per pack:

Disc
 J75-1450
 RA

Cone
 J75-1452
 RA

Disc
 J75-1454
 RA

Cone
 J75-1456
 RA

\$30.00

Acrylic Polishers

Three-step polishing system for removing, finishing and high luster polishing of acrylics. Silicon carbide embedded in a medium-soft, flexible synthetic rubber matrix.

6 per Pkg.
 Rec. RPM:
 5,000-7,000



Desc:
 Shank:

Small Point
 HP

Point
 HP

Large Point
 HP

Jumbo
 HP

Shaping / Coarse

Item Code:

J75-1392

J75-1394

J75-1396

J75-1398

Smoothing / Medium

Item Code:

J75-1400

J75-1402

J75-1404

J75-1406

High Shine / Fine

Item Code:

J75-1408

J75-1410

J75-1412

J75-1414

Price:

\$17.25

Green & White Stones

Green Stones contain silicon carbide and are excellent for fast contouring and finishing of porcelain, composites, and gold and silver alloys. White Stones, made of aluminum oxide, are ideal for fine contouring and finishing of enamel, composites, and porcelains.

12 per Pkg.
 Rec. RPM:
 8,000-12,000



Desc:
 Item Code:
 Shank:
 Price:

WH-6
 J75-1416
 HP

TC-4
 J75-1418
 HP

FL-3
 J75-1420
 HP

IC-4
 J75-1422
 HP

CN-1
 J75-1421
 FG

FL-2
 J75-1423
 FG

\$13.95



Desc:

IC-3

IC-5

IC-7

IC-9

KN-7

PC-2

WH-3

WH-6

CY-3

TC-2

TC-4

CN-1

FL-2

Item Code:

J75-1430

J75-1432

J75-1434

J75-1436

J75-1438

J75-1440

J75-1442

J75-1444

J75-1424

J75-1426

J75-1428

J75-1445

J75-1446

Shank:

HP

HP

HP

HP

HP

HP

HP

HP

HP

HP

HP

FG

FG

Price:

\$13.95