

In-office digital impressioning & milling solutions

In-office digital impression & milling package value package **\$67,995⁰⁰** (1-2-3-4-5)

Package includes:

1 **MEDIT**
i-500 Intra
Oral Scanner



High ROI

The design and soft ware allow for an unparalleled performance with a very competitive price.

Vast Flexibility

The i500 uses an open CAD/CAM system, allowing you to export STL files.

Impressive Speed

Using two high-speed cameras, the i500 scans quickly and efficiently.

Powderless

The i500 has no need for the use of powder for regular cases, making the experience more convenient and comfortable for the patient.

Small Tip

The small tip on the i500 is beneficial for both the user and the patient through simpler maneuver ability and improved patient comfort.

Video Type

With the i500, there is no need to worry about choppy results, as it takes rapid video - based scans.

Color

The i500 scanner produces precise-colored scans which can distinguish between soft tissue, plaque and teeth.

3 **exocad** | chairside CAD software

More automation, simplified user experience, Optimized for usage in a clinical environment. Includes Lap top PC with package.

2 **vhf** | Z4 milling machine

- Most finely balanced high frequency spindle with up to 100,000 RPM
- Precise guiding rails and high-class ballscrew drives
- Massive internal machine bed made of aluminium cast
- Premium industrial quality
- Uses a large variety of materials: Glass Ceramics, PMMA, Zirconium Oxide, Composites, Pre fabricated Titanium Abutments



4 Package training

Receive FREE

5 **MoonRAY**
3D printing machine

Your choice:
model S100 or D75

Retail value:
\$3995



D75

S100