




# Conch series

Rate:30V/DC  
 Frequency:28KHz±3KHz  
 ENDO function is optional



# Coral series

Rate:30V/DC  
 Frequency:28KHz±3KHz

 <b>P5L</b> L3 Handpiece SCALING / PERIO	 <b>P5</b> H1 Handpiece SCALING / PERIO	 <b>BD2</b> H3 Handpiece SCALING / PERIO
 <b>P4</b> H2 Handpiece SCALING / PERIO	 <b>B5L</b> L1 Handpiece SCALING / PERIO	 <b>B5/B55</b> H1 / H3 Handpiece SCALING / PERIO

 <b>P7L</b> L3 Handpiece SCALING / PERIO / ENDO	 <b>P7</b> H1 Handpiece SCALING / PERIO / ENDO	 <b>P6</b> H2 Handpiece SCALING / PERIO / ENDO
--	---	---