

Pearl series P9L

Auto-water Supply Scaler

SCALING/ ENDO/ PERIO

- Auto-water supply device;
- Built-in power supply;
- Detachable handpiece can be autoclaved under the high temperature and pressure (135°C, 0.22MPa)
- Aluminum alloy LED handpiece;
- Obtaining national patents;
- Identifying scaler tip automatically;
- PC control, automatic frequency tracking.



Compact, space saving: 43x145x150(mm)

Waterway system without consumables, cost saving.

White light is optional



Detachable aluminum alloy LED handpiece.

Noiseless working

No Consumables!



P9 CE

Handpiece: H3



P9-2

Handpiece: D2



P9-3

Handpiece: D3