






CARMEL

www.carmelindustries.com

DENTAL WAX CATALOG





Bite Block
Packaged 120 units per box.

BB40019 Pink




Bite Rim
Packaged 100 units per box.

BB40219 Pink




Occlusal Rim
Packaged 100 units per box.

BB40319 Pink



Bite Stick
2 lengths available: 9" (225 mm) or 4 1/2" (115 mm), width and thickness 3/8" (10 mm). Available in soft, medium and hard. Available in pink or yellow.
Packaged 4 lb (1.85 kg) per box.

4.5" x 3/8" x 3/8" - 115 mm x 10 mm x 10 mm
BS40019 Pink
BS40002 Yellow
9" x 3/8" x 3/8" - 225 mm x 10 mm x 10 mm
BS40119 Pink
BS40102 Yellow



5 LBS	1 LB	5 LBS	1 LB	5 LBS	1 LB	5 LBS	1 LB	5 LBS	1 LB
Universal		Q-Wax		100 (Regular)		101 (Thin)		102 (Thick)	
BP40009	BP40109	BP40409	BP40509	BP44409	BP44509	BP44609	BP44809	BP44609	BP44619
BP40019	BP40119	BP40419	BP40519	BP44419	BP44519	BP44619	BP44619	BP44619	BP44619
Extra Tough		Bulldog (Premium Tenacious)		TruForm		Flexi-Tough		Setup Wax	
BP42519	BP42619	BP47319	BP47419	BP42109	BP42209	BP43019	BP43119	BP40219	BP40319
BP42119	BP42219								
All Season		Ultra Precision		Bite Wax		Pink #2		Pink #3	
BP44009	BP44109	BP45619	BP45719	BP43402	BP43502	BP40619	BP40719	BP40819	BP40919
BP44019	BP44119								
Pretty Pink Soft		Pretty Pink Medium		Pretty Pink Hard		Pretty Pink X-Hard		Master Wax	
BP41219	BP41319	BP41419	BP41519	BP41619	BP41717	BP41819	BP41919		BP42009
Beeswax Sheets									
BP43290	BP43390								

Baseplate
(Dimensional)
throughout the American continent
All formulas are available in 1 lb (454 g)
and 5 lb (2.270 kg) boxes.
5 lb boxes packaged 10 boxes/case

Sheet dimensions
3 inch - 7.62 cm
6 inch - 15.24 cm

Pink
Red
Yellow



Modelling Wax Sheets (European dimensions & formulas)

All formulas are available in 1.1 lb (500 g) and 5.51 lb (2.5 kg) boxes.
2.5 kg boxes packaged 10 boxes/case.

Sheet dimensions
3 5/16 inch - 8.4 cm
6 3/4 inch - 17.14 cm

2.5 KG	500 G	2.5 KG	500 G	2.5 KG	500 G
Ultra Precision		Summit (Regular)		Tough #4	
ME40619	ME40719	ME40219	ME40319	ME40019	ME40119
Setup Wax Soft		Setup Wax Regular		Setup Wax Hard	
ME43019	ME44019	ME45019	ME46019	ME47019	ME48019
Euro #2		Summit (Hard)			
	ME40919	ME40419	ME40519		

Pink



Presentation Wax
Packaged in 2.46 oz (70 g) tins.

 PW40003 White




Chrome Cobalt Wax
1 lb (454 g) of chunks per box.

 IW42014 Dark Blue






Sculpturing Wax
Packaged in 2.46 oz (70 g) tins.

-  SC40001 Purple
-  SC40004 Blue
-  SC40005 Brown
-  SC40009 Red
-  SC40010 Green
-  SC40018 Gray







Modelling Wax
Packaged in 2.46 oz (70 g) tins.

-  IW42104 Blue
-  IW42109 Red
-  IW42110 Green



Margin Wax
Packaged in 2.46 oz (70 g) tins.

-  MW40002 Yellow
-  MW40003 White
-  MW40009 Red
-  MW40014 Dark Blue



Inlay Wax

Tins: 2.46 oz (70 g) per tin.
 Chunks: 1 lb (454 g) per box,
 packaged 24 boxes/case.
 Sticks: 12, 24, or 120 per box.


TINS				
	IW41004	Blue		
	IW41010	Green		
	IW41017	Ivory		
CHUNKS				
	IW40004	Blue		
	IW40010	Green		
	IW40017	Ivory		
STICKS				
	12x	24x	120x	
	IW40104	IW40204	IW40304	Blue
	IW40110	IW40210	IW40310	Green
	IW40117	IW40217	IW40317	Ivory



Sticky Wax

Chunks: 1 lb (454 g) per box, packaged 24 boxes/case.
 Sticks: 12, 24, or 120 per box.

CHUNKS HONEY			
	SW40007		
STICKS HONEY			
	12x	24x	120x
	SW40207	SW40307	SW40407



0.5 LB. SPOOLS				250 G SPOOLS		
	Blue				Blue	
Gauge	Round	1/2 Round	1/2 Pear	Metric	Round	1/2 Round
4 (0,204", 5.2mm)	SP42004	SP41004	SP41704	0.5 mm	SP43004	
6 (0,162", 4.1mm)	SP40004	SP41104		1.0 mm	SP43104	
8 (0,128", 3.3mm)	SP40104	SP41204		1.5 mm	SP43204	SP44204
10 (0,102", 2.6mm)	SP40204	SP41304		2.0 mm	SP43304	SP44304
12 (0,081", 2.1mm)	SP40304	SP41404		2.5 mm	SP43404	SP44404
14 (0,064", 1.6mm)	SP40404	SP41504		3.0 mm	SP43504	SP44504
16 (0,051", 1.3mm)	SP40504	SP41604		3.5 mm	SP43604	SP44604
18 (0,040", 1.0mm)	SP40604			4.0 mm	SP43704	
20 (0,032", 0.8mm)	SP40704			4.5 mm	SP43804	
22 (0,025", 0.6mm)	SP40804			5.0 mm	SP43904	
24 (0,020", 0.5mm)	SP40904			5.5 mm	SP44004	
				6.0 mm	SP44104	

Wire Wax

Wax wire wound on spools containing 0.5 lb (227 g) or 8.8 oz (250 g).

Available in blue, natural, green, turquoise and white, in different hardnesses.

Packaged 1 spool/box, 50 boxes/case.
 Different sizes can be combined to fill a case.



Dipping Wax

Packaged in 9 oz (255 g) and 18 oz (510 g) bottles.

9 OZ		18 OZ	
	DW40002	DW40102	Yellow
	DW40004	DW40104	Blue
	DW40009	DW40109	Red
	DW40010	DW40110	Green



Utility Wax


Two profiles available:

Round: 3/16" (4.76 mm) diameter x 11" (27.94 cm) long. 64 rods per box.

Square: 3/16" (4.76 mm) x 11" (27.94 cm) long.

44 strips per box, packaged 50 boxes/case.

ROUND	SQUARE	
	UW40103	White
	UW40109	Red



Boxing Wax

11" (27.94 cm) long x 1.5" (3.81 cm) wide

Regular: 0.062" (1.6 mm) thick.


Extra thin: 0.040" (1.0 mm) thick.

Packaged 1 lb (454 g) & 5 lb (2.270 kg) per box.

1 lb boxes packaged 50 boxes/case.

5 lb boxes packaged 10 boxes/case.

	REGULAR	EXTRA THIN	
	1 lb BW40309	BW40409	Red
	5 lbs BW40509	BW40609	Red



Periphery Wax

6" (15.24 cm) x 1/4" (6.4 mm) diameter.

Packaged 60 rods per box.

PW40004 Blue



Lingual Bar

Packaged 1 lb (454g) per box.

20 boxes/case.

LW40019 Pink




Palate Wax Forms

Prefabricated wax patterns.

Packaged 80 units/box.

LW40019 Pink



Sprue Wax

6" (15.24 cm) long round rod.

Packaged 1 lb (454 g) per box, 20 boxes/case.

Diameters | Diamètres | Diámetros

5 gauge (0.182", 4.62mm), 6 gauge (0.162", 4.1mm), 8 gauge (0.128", 3.3mm),

10 gauge (0.102", 2.6mm) and 1/4" inch (6.4mm).

	GAUGE				INCH	
	5	6	8	10	1/4	Red
	SW40009	SW40109	SW40209	SW40309	SW40409	



Beading Wax

Packaged 5 lb (2.27 kg) per box,

10 boxes/case.

BW41009 Red



Non-Laminated Bite Wafers

Single layered bite wafer in U form.
Packaged 50 units per box.

Laminated Bite Wafers

Double layered bite wafer in U form stabilized by an aluminium foil in the center.
Packaged 24, 36 or 84 units/box, 16 boxes/case.

NON-LAMINATED BITE WAFER

Soft 50/box

Hard 50/box

 BW40002
Yellow

 BW40010
Green

LAMINATED BITE WAFER

24/box

36/box

84/box


 BW40004
Blue

 BW40005
Cooper

 BW40105
Cooper


 BW40003
White

 BW40103
White



Orthopedic Wax Bites

Packaged 28 units per box, 21 boxes/case.



Bite Registration Wax

2 forms available:

Sheets: Dimensions: 6" (15.24 cm) x 3" (7.62 cm).

Packaged 1 lb (454 g)/box


Wafers: Packaged 25 per box).

 OW40019 Pink

SHEETS | WAFERS

 BR40004 BR40104 Blue

Used by orthodontic patients to guard against irritation.



Relief Wax
(Unscented box containing plain wax)

An unscented polypropylene box containing plain wax.

Packaged: 50 Relief Wax plastic boxes/box, 40 boxes/case; or 100 Relief Wax plastic boxes/box, 20 boxes/case.

-  RW40180 Translucent box
-  RW40113 Silver box



Relief Wax (scented colored box)

Scented colored polypropylene boxes containing plain wax. Packaged: 50 Relief Wax plastic boxes/box, 40 boxes/case or 100 Relief Wax plastic boxes/box, 20 boxes/case.

Scented colored box

-  RW40201 Grape
-  RW40202 Pineapple
-  RW40203 Coconut
-  RW40204 Blueberry
-  RW40205 Chocolate
-  RW40209 Cherry
-  RW40210 Mint
-  RW40212 Orange
-  RW40216 Cinnamon
-  RW40219 Bubble Gum
-  RW40221 Cotton Candy



Orthodontic Wax (transparent box)

A transparent box made of crystal clear polystyrene. Available with plain or mint wax. 100 Orthodontic Wax plastic boxes/box.

-  RW40080 Plain wax
-  RW40010 Mint wax



Orthodontic Wax Blister Pack

A blister pack containing a single Orthodontic wax box (a transparent box made of crystal clear polystyrene). Two distinct blister packs available for Orthodontic wax boxes containing plain or mint wax.

-  RW40480 Plain wax
-  RW40410 Mint wax



Patient Wax (flavored colored wax)

A translucent polypropylene box containing flavored colored wax.


Packaged: 50 Relief Wax plastic boxes/box, 40 boxes/case; or 100 Relief Wax plastic boxes/box, 20 boxes/case.

Flavored colored wax

-  RW40312 Orange
-  RW40323 Lime
-  RW40301 Grape
-  RW40302 Lemon
-  RW40309 Cherry

CARMEL

CARMEL'S OTHER WAX PRODUCT CATALOG



GENERAL WAXES CATALOG



WWW.CARMELINDUSTRIES.COM



HEAD OFFICE

CARMEL GROUP INC.

10220 Armand-Lavergne Ave.
Montreal, Quebec, Canada H1H 3N5
Tel. 514 270-5377 | Fax. 514 270-2025

info@carmelindustries.com
www.carmelindustries.com

U.S. ADDRESS

CARMEL GROUP INC.

P.O. Box 3150
Champlain, New York
USA 12919

CARMEL