

Color numbering and coding



	Grain	Grain Size	ISO Indication	
	UF	Ultra Fine	4-14 μ	494
	C	Super Fine	10-36 μ	504
	F	Fine	27-76 μ	514
	/	Medium - Standard	66-126 μ	524
	G	Coarse	108-181 μ	534
	SG	Super Coarse	151-213 μ	544

Shank sizes and coding

\emptyset		Length	ISO Indication	
1.6	FGS	Friction Grip Short	16mm	313
1.6	FG	Friction Grip	19mm	314
1.6	FGXL	Friction Grip Long	25mm	316
1.6	FGXXL	Friction Grip Extra Long	30mm	317
2.35	HP	Hand Piece	44.5mm	104
2.35	HPL	Hand Piece Long	65mm	105
2.35	HPXL	Hand Piece Extra Long	70mm	106
2.35	RA	Right Angle	22mm	204
2.35	RAL	Right Angle Long	26mm	205
2.35	RAXL	Right Angle Extra Long	34mm	206

Packing information

S/C.U-13.020.RAXL-5	"S" indicates Sterile
S/C.FD6.010.FG	"S" indicates Sterile
S.1.018.RA	"S" here is indicating Steel!
P.ZR5.130.HP-1	130 is for radius and "-1" is for packing
P.U8M.RA-10	"-10" indicates 10 piece pack
P.RM734.090.HP/E-12	"-12" is for 12 piece packing
P.PRONYL5-PINK.RA-100	"-100" is for 100 piece packing
P.FLEX/POPON.G.140	"140" is for radius
P.FLEX/POPON.M.100	"100" is for radius
M.4011.HP/E-6	"-6" is for 6 piece packing

Frank Dental product code indications

D	Diamond
C	Carbide
S	Steel
P	Polisher
M	Mandrel
ZR	Zircon
S/C	Sterile Carbide
S/D	Sterile Diamond