

E-Sapphire Ceramic Brackets



>>>>>

- Transparent aesthetic bracket, polish bracket surface for three times elegant and generous
- Low wing, low friction, approximate metal bracket
- CNC technology, more smooth wing edge, more perfect radian and comfortable wearing
- Bottom laser maintenance bead greatly improve bonding strength
- Good concealing, uniformly forced, crystal clear, especially suitable for metal allergy patient



E-Sapphire Ceramic Bracket-Roth

Item No. : ES-101

	Bracket	Torq.	Ang.	Width(mm)	0.018R	0.018L	0.022R	0.022L
Maxillary	Central	+12°	+ 5°	3.60	ES-101-811	ES-101-821	ES-101-211	ES-101-221
	Lateral	+ 8°	+ 9°	3.10	ES-101-812	ES-101-822	ES-101-212	ES-101-222
	Cuspid	- 2°	+13°	3.20	ES-101-813	ES-101-823	ES-101-213	ES-101-223
	Cuspid W/H	- 2°	+13°	3.20	ES-101-813H	ES-101-823H	ES-101-213H	ES-101-223H
	1st Bicuspид	- 7°	0°	3.20	ES-101-814	ES-101-824	ES-101-214	ES-101-224
	1st Bicuspид W/H	- 7°	0°	3.20	ES-101-814H	ES-101-824H	ES-101-214H	ES-101-224H
	2nd Bicuspид	- 7°	0°	3.20	ES-101-815	ES-101-825	ES-101-215	ES-101-225
	2nd Bicuspид W/H	- 7°	0°	3.20	ES-101-815H	ES-101-825H	ES-101-215H	ES-101-225H
Mandibular	Central	- 1°	0°	2.60	ES-101-841	ES-101-831	ES-101-241	ES-101-231
	Lateral	- 1°	0°	2.60	ES-101-842	ES-101-832	ES-101-242	ES-101-232
	Cuspid	-11°	+7°	3.20	ES-101-843	ES-101-833	ES-101-243	ES-101-233
	Cuspid W/H	-11°	+7°	3.20	ES-101-843H	ES-101-833H	ES-101-243H	ES-101-233H
	1st Bicuspид	-17°	0°	3.20	ES-101-844	ES-101-834	ES-101-244	ES-101-234
	1st Bicuspид W/H	-17°	0°	3.20	ES-101-844H	ES-101-834H	ES-101-244H	ES-101-234H
	2nd Bicuspид	-22°	0°	3.20	ES-101-845	ES-101-835	ES-101-245	ES-101-235
	2nd Bicuspид W/H	-22°	0°	3.20	ES-101-845H	ES-101-835H	ES-101-245H	ES-101-235H

Slot	Full Set Item Number		
0.018	ES-101-8A (3.4.5# with hook)	ES-101-8B (3# with hook)	ES-101-8C (no hook)
0.022	ES-101-2A (3.4.5# with hook)	ES-101-2B (3# with hook)	ES-101-2C (no hook)

E-Sapphire Ceramic Bracket-MBT

Item No. : ES-102

	Bracket	Torq.	Ang.	Width(mm)	0.018R	0.018L	0.022R	0.022L
Maxillary	Central	+17°	+4°	3.60	ES-102-811	ES-102-821	ES-102-211	ES-102-221
	Lateral	+10°	+8°	3.10	ES-102-812	ES-102-822	ES-102-212	ES-102-222
	Cuspid	0°	+8°	3.20	ES-102-813	ES-102-823	ES-102-213	ES-102-223
	Cuspid W/H	0°	+8°	3.20	ES-1012-813H	ES-102-823H	ES-102-213H	ES-102-223H
	1st Bicuspид	- 7°	0°	3.20	ES-102-814	ES-102-824	ES-102-214	ES-102-224
	1st Bicuspид W/H	- 7°	0°	3.20	ES-102-814H	ES-102-824H	ES-102-214H	ES-102-224H
	2nd Bicuspид	- 7°	0°	3.20	ES-102-815	ES-102-825	ES-102-215	ES-102-225
	2nd Bicuspид W/H	- 7°	0°	3.20	ES-102-815H	ES-102-825H	ES-102-215H	ES-102-225H
Mandibular	Central	- 6°	0°	2.60	ES-102-841	ES-102-831	ES-102-241	ES-102-231
	Lateral	- 6°	0°	2.60	ES-102-842	ES-102-832	ES-102-242	ES-102-232
	Cuspid	0°	+3°	3.20	ES-101-843	ES-102-833	ES-102-243	ES-102-233
	Cuspid W/H	0°	+3°	3.20	ES-102-843H	ES-102-833H	ES-102-243H	ES-102-233H
	1st Bicuspид	-12°	+2°	3.20	ES-102-844	ES-102-834	ES-102-244	ES-102-234
	1st Bicuspид W/H	-12°	+2°	3.20	ES-102-844H	ES-102-834H	ES-102-244H	ES-102-234H
	2nd Bicuspид	-17°	+2°	3.20	ES-102-845	ES-102-835	ES-102-245	ES-102-235
	2nd Bicuspид W/H	-17°	+2°	3.20	ES-102-845H	ES-102-835H	ES-102-245H	ES-102-235H

Slot	Full Set Item Number		
0.018	ES-102-8A (3.4.5# with hook)	ES-102-8B (3# with hook)	ES-102-8C (no hook)
0.022	ES-102-2A (3.4.5# with hook)	ES-102-2B (3# with hook)	ES-102-2C (no hook)

Metal Self-ligating Brackets

Ceramic Self-ligating Brackets

Metal Brackets

Ceramic Brackets

Buccal Tubes

Molar Bands

Attachments & Archwires

Elastomerics

Instruments