


E-Unique Brackets



>>>>>
Low profile means more comfortable

25% Less



E-Unique



Common Style

E-Unique Bracket-Roth Mini Bondable

Item No. : EM-301

	Bracket	Torq.	Ang.	Width(mm)	0.018R	0.018L	0.022R	0.022L
Maxillary	Central	+12°	+ 5°	3.58	EM-301-811	EM-301-821	EM-301-211	EM-301-221
	Lateral	+ 8°	+ 9°	2.95	EM-301-812	EM-301-822	EM-301-212	EM-301-222
	Cuspid	- 2°	+10°	3.06	EM-301-813	EM-301-823	EM-301-213	EM-301-223
	Cuspid W/H	- 2°	+10°	3.06	EM-301-813H	EM-301-823H	EM-301-213H	EM-301-223H
	1st Bicuspid	- 7°	0°	3.08	EM-301-814	EM-301-824	EM-301-214	EM-301-224
	1st Bicuspid W/H	- 7°	0°	3.08	EM-301-814H	EM-301-824H	EM-301-214H	EM-301-224H
	2nd Bicuspid	- 7°	0°	3.08	EM-301-815	EM-301-825	EM-301-215	EM-301-225
	2nd Bicuspid W/H	- 7°	0°	3.08	EM-301-815H	EM-301-825H	EM-301-215H	EM-301-225H
Mandibular	Central	- 1°	0°	2.55	EM-301-841	EM-301-831	EM-301-241	EM-301-231
	Lateral	- 1°	0°	2.55	EM-301-842	EM-301-832	EM-301-242	EM-301-232
	Cuspid	-11°	+7°	3.12	EM-301-843	EM-301-833	EM-301-243	EM-301-233
	Cuspid W/H	-11°	+7°	3.12	EM-301-843H	EM-301-833H	EM-301-243H	EM-301-233H
	1st Bicuspid	-17°	0°	3.12	EM-301-844	EM-301-834	EM-301-244	EM-301-234
	1st Bicuspid W/H	-17°	0°	3.12	EM-301-844H	EM-301-834H	EM-301-244H	EM-301-234H
	2nd Bicuspid	-22°	0°	3.12	EM-301-845	EM-301-835	EM-301-245	EM-301-235
	2nd Bicuspid W/H	-22°	0°	3.12	EM-301-845H	EM-301-835H	EM-301-245H	EM-301-235H

Slot	Full Set Item Number			
0.018	EM-301-8A (3.4.5# with hook)	EM-301-8B (3# with hook)	EM-301-8C (no hook)	
0.022	EM-301-2A (3.4.5# with hook)	EM-301-2B (3# with hook)	EM-301-2C (no hook)	

E-Unique Bracket-MBT Mini Bondable

Item No. : EM-302

	Bracket	Torq.	Ang.	Width(mm)	0.018R	0.018L	0.022R	0.022L
Maxillary	Central	+17°	+4°	3.55	-	-	EM-302-211	EM-302-221
	Lateral	+10°	+8°	2.93	-	-	EM-302-212	EM-302-222
	Cuspid	- 7°	+8°	3.10	-	-	EM-302-213	EM-302-223
	Cuspid W/H	- 7°	+8°	3.10	-	-	EM-302-213H	EM-302-223H
	1st Bicuspid	- 7°	0°	3.05	-	-	EM-302-214	EM-302-224
	1st Bicuspid W/H	- 7°	0°	3.05	-	-	EM-302-214H	EM-302-224H
	2nd Bicuspid	- 7°	0°	3.05	-	-	EM-302-215	EM-302-225
	2nd Bicuspid W/H	- 7°	0°	3.05	-	-	EM-302-215H	EM-302-225H
Mandibular	Central	- 6°	0°	2.55	-	-	EM-302-241	EM-302-231
	Lateral	- 6°	0°	2.55	-	-	EM-302-242	EM-302-232
	Cuspid	- 6°	+3°	3.10	-	-	EM-302-243	EM-302-233
	Cuspid W/H	- 6°	+3°	3.10	-	-	EM-302-243H	EM-302-233H
	1st Bicuspid	-12°	+2°	3.05	-	-	EM-302-244	EM-302-234
	1st Bicuspid W/H	-12°	+2°	3.05	-	-	EM-302-244H	EM-302-234H
	2nd Bicuspid	-17°	+2°	3.05	-	-	EM-302-245	EM-302-235
	2nd Bicuspid W/H	-17°	+2°	3.05	-	-	EM-302-245H	EM-302-235H

Slot	Full Set Item Number			
0.018	-	-	-	
0.022	EM-302-2A (3.4.5# with hook)	EM-302-2B (3# with hook)	EM-302-2C (no hook)	

Metal Self-ligating Brackets

Ceramic Self-ligating Brackets

Metal Brackets

Ceramic Brackets

Buccal Tubes

Molar Bands

Attachments & Archwires

Elastomerics

Instruments