


COMPOSITE MICROFILL POLISHERS

# KENDA COMPOSITE MICROFILL


## 2-step polishing system

Top selling polisher for all microfilled composite restorations. Achieves an optimal and long lasting high shine.

Also suited for simply removing cement remnants of orthodontic bands (brackets).



Set REF 750



<b>REF</b> reorder	4001	4002	4003	4004	4005	4006	4007	4008	5001	5002	5003	5004	5005	5006	5007	5008	
<b>ISO</b> 638 204...	243 523 055	254 523 050	304 523 100	030 523 050	030 523 060	243 523 030	237 523 035	243 523 050	638 204...	243 493 055	254 493 050	304 493 100	030 493 050	030 493 060	243 493 030	237 493 035	243 493 050
<b>ISO</b> Size Ø×L mm	5,4×13	5×8	10×2	5×8	6,5×10	3×7	3×6	5×10	5,4×13	5×8	10×2	5×8	6,5×10	3×7	3×6	5×10	
<b>μ</b> Grit	medium								ultra-fine								
<b>Color</b>	yellow								blue								
<b>rec.</b>	7 500 - 10 000 RPM								5 000 - 7 500 RPM								
<b>max.</b>	10 000 RPM								10 000 RPM								

- Box of: 6, 12, 25, 100
- Blister of: 6×1, 6×2, 72×1, 72×2
- N** Pressure: 0,3 - 2 N
- Autoclavable: 135 °C, 275 °F
- Shank: RA, stainless steel
- Medical device: class IIa