



## Membrane Fixed Screw Set (Screw 54EA)

Consist of short length screws(3 & 4 mm) for membrane fixation(Ti-mesh).



### Fixing Membrane Screw

Code	D*L	Quantity
AU-12-003	Ø1.2 X 3	18EA
AU-12-004	Ø1.2 X 4	18EA
AU-14-004	Ø1.4 X 4	18EA



Membrane Fixed Screw Set


TMS-02



use for handpiece use for handle


## Block Bone Screw Set (Screw 54EA)

Consist of long length screws (8, 10, 12mm) for block bone fixation or Tenting purpose.



### Block Bone Screw

Code	D*L	Quantity
AU-16-008	Ø1.6 X 6	6EA
AU-16-010	Ø1.6 X 8	6EA
AU-16-012	Ø1.6 X 10	6EA
Code	D*L	Quantity
AU-18-008	Ø1.8 X 6	6EA
AU-18-010	Ø1.8 X 8	6EA
AU-18-012	Ø1.8 X 10	6EA
Code	D*L	Quantity
AU-20-008	Ø2.0 X 6	6EA
AU-20-010	Ø2.0 X 8	6EA
AU-20-012	Ø2.0 X 10	6EA



Block Bone Screw Set

TMS-03



use for handpiece use for handle

## Size of Bone Screw

Code	D*L	Code	D*L	Code	D*L	Code	D*L	Code	D*L
AU-12-003	Ø1.2 X 3								
AU-12-004	Ø1.2 X 4	AU-14-004	Ø1.4 X 4	AU-16-004	Ø1.6 X 4	AU-18-004	Ø1.8 X 4	AU-20-004	Ø2.0 X 4
AU-12-006	Ø1.2 X 6	AU-14-006	Ø1.4 X 6	AU-16-006	Ø1.6 X 6	AU-18-006	Ø1.8 X 6	AU-20-006	Ø2.0 X 6
AU-12-007	Ø1.2 X 7	AU-14-007	Ø1.4 X 7	AU-16-007	Ø1.6 X 7	AU-18-007	Ø1.8 X 7	AU-20-007	Ø2.0 X 7
AU-12-018	Ø1.2 X 8	AU-14-008	Ø1.4 X 8	AU-16-008	Ø1.6 X 8	AU-18-008	Ø1.8 X 8	AU-20-008	Ø2.0 X 8
AU-12-010	Ø1.2 X 10	AU-14-010	Ø1.4 X 10	AU-16-010	Ø1.6 X 10	AU-18-010	Ø1.8 X 10	AU-20-010	Ø2.0 X 10
AU-12-012	Ø1.2 X 12	AU-14-012	Ø1.4 X 12	AU-16-012	Ø1.6 X 12	AU-18-012	Ø1.8 X 12	AU-20-012	Ø2.0 X 12