


# Piezo Ultrasonic Scaler Units

## Woodpecker

### Piezo Ultrasonic Scaler Replacement Tips

#### Periodontics




*for subgingival and supragingival calculus.*

**Stainless Steel**

**P1 (EMS® type)**.....G09-0030

**PD1 (Satelec® type)**.....G09-0060



*for fine scaling and deep perio calculus.*

**Stainless Steel**

**P3 (EMS® type)**.....G09-0032

**PD3 (Satelec® type)**.....G09-0062




*for smoothing overhangs and debride in-depth root (diamond-coated)*

**Stainless Steel**

**P3D (EMS® type)**.....G09-0033

**PD3D (Satelec® type)**.....G09-0063



*for subgingival and interdental calculus.*

**Stainless Steel**

**P4 (EMS® type)**.....G09-0034

**PD4 (Satelec® type)**.....G09-0064




*for periodontal scaling*

**Stainless Steel**

**P2L (EMS® type)**.....G09-0116

**PD2L (Satelec® type)**.....G09-0282



*for periodontal scaling*

**Stainless Steel**

**P2R (EMS® type)**.....G09-0118

**PD2R (Satelec® type)**.....G09-0286

#### Scaling



*for supragingival calculus and plaque. Neck of teeth and hard calculus on the teeth point.*

**Stainless Steel**

**G1 (EMS® type)**.....G09-0016

**GD1 (Satelec® type)**.....G09-0046

**Titanium**

**G1T (for EMS® type)**.....G09-0160



*for supragingival heavy calculus.*

**Stainless Steel**

**G2 (EMS® type)**.....G09-0018

**GD2 (Satelec® type)**.....G09-0048

**Titanium**

**G2T (for EMS® type)**.....G09-0162



*for subgingival, supragingival, interdental calculus and plaque.*

**Stainless Steel**

**G3 (EMS® type)**.....G09-0020

**GD3 (Satelec® type)**.....G09-0050



*for supragingival calculus and plaque.*

**Stainless Steel**

**G4 (EMS® type)**.....G09-0022

**GD4 (Satelec® type)**.....G09-0052

**Titanium**

**GD4T (Satelec® type)**.....G09-0444



*for supragingival, neck calculus and plaque*

**Stainless Steel**

**G5 (EMS® type)**.....G09-0024

**GD5 (Satelec® type)**.....G09-0054

**Titanium**

**GD5T (Satelec® type)**.....G09-0446



*for supragingival heavy calculus*

**Stainless Steel**

**G6 (EMS® type)**.....G09-0026

**GD6 (Satelec type)**.....G09-0056

#### Cavity Preparation



*For occlusal surface and cervical margins preparation.*

**Stainless Steel**

**SB1 (EMS® type)**.....G09-0124

**SBD1 (Satelec® type)**.....G09-0296



*Remove the carious teeth without damaging the adjacent teeth*

**Stainless Steel**

**SB2 (EMS® type)**.....G09-0126

**SBD2 (Satelec® type)**.....G09-0298



*Remove the carious teeth without damaging the adjacent teeth*

**Stainless Steel**

**SB3 (EMS® type)**.....G09-0128

**SBD3 (Satelec® type)**.....G09-0300