



ROUND




1/4 1/2 1 2 3 4 5 6 7 8 10



INVERTED CONE



33 1/2 34 35 36 37 37L 38 39



PEAR



1/2 P 1P 2P 4P 329 330 300M 331 332 330L 331L 332L

FISSURE- FLAT END CROSS CUT




556 557 558 559 560 556L 557L 558L

FISSURE- ROUND END CROSSCUT



1556 1557 1558

FISSURE- FLAT END PLAIN




56 57 57L 58 59

FISSURE- ROUND END PLAIN



1156 1157 1158

TAPER/FLAT END PLAIN



169 170 171 169L 170L 171L

AMALGAM PREP



244 245 246 256 257 271

TAPER/ FLAT END CROSS CUT



699 700 701 702 703 699L 700L 701L 702L

TAPER/ ROUND END PLAIN



1170 1171 1170L 1171L

END CUT



956 957

WHEEL



14