


**W368**  
Pointed Bulb  
3




Shank	L mm	ISO	Order No.		
FG		• W806 314 257 524...	W 368.314...	020	
		• W806 314 257 514...	FW 368.314...	020	
		• W806 314 257 534...	GW 368.314...	020	023

**W379**  
Egg Shape / Football  
3




Shank	L mm	ISO	Order No.		
FG		• W806 314 277 524...	W 379.314...	023	
		• W806 314 277 514...	FW 379.314...	023	
		• W806 314 277 534...	GW 379.314...	023	

**W801**  
Round  
3




Shank	L mm	ISO	Order No.			
FG		• W806 314 001 524...	W 801.314...	014	021	023
		• W806 314 001 514...	FW 801.314...	014	021	
		• W806 314 001 534...	GW 801.314...	014	021	023

**W811**  
Barrel  
3




Shank	L mm	ISO	Order No.		
FG		• W806 314 038 534...	GW 811.314...	033	

**W830L**  
Pear long  
3




Shank	L mm	ISO	Order No.			
FG		• W806 314 239 524...	W 830L.314...	012	014	016
		• W806 314 239 514...	FW 830L.314...	012	014	016
		• W806 314 239 534...	GW 830L.314...	012	014	016

**W835KR**  
Round Edge Cylinder  
3



Shank	L mm	ISO	Order No.		
FG		• W806 314 156 524...	W 835KR.314...	012	
		• W806 314 156 514...	FW 835KR.314...	012	
		• W806 314 156 534...	GW 835KR.314...	012	

**W837KR**  
Round Edge Cylinder  
3




Shank	L mm	ISO	Order No.		
FG		• W806 314 158 524...	W 837KR.314...	014	
		• W806 314 158 514...	FW 837KR.314...	014	
		• W806 314 158 534...	GW 837KR.314...	014	
		• W806 314 158 544...	SGW 837KR.314...		016

**W845KR**  
Round Edge Taper  
3




Shank	L mm	ISO	Order No.		
FG		• W806 314 544 524...	W 845KR.314...	016	

**W846KR**  
Round Edge Taper  
3







Shank	L mm	ISO	Order No.		
FG		• W806 314 545 524...	W 846KR.314...	016	023

**W847**  
Flat End Taper  
3



Shank	L mm	ISO	Order No.		
FG		• W806 314 172 524...	W 847.314...		016
		• W806 314 172 534...	GW 847.314...	014	016

 Cavity preparation	 Crown preparation	 Filling process	 Jaw surgery
--	---	---	--

- blue ring standard – 126µm, ISO No. 524
- red ring fine – 40µm, ISO No. 514
- green ring coarse – 151µm, ISO No. 534
- black super-coarse – 181µm, ISO No. 544

**W847KR**

Round End Taper

3



Shank	L mm	ISO	Order No.	8,0	8,0
FG		• W806 314 546 524...	W 847KR.314...	016	
		• W806 314 546 514...	FW 847KR.314...	016	
		• W806 314 546 534...	GW 847KR.314...	016	018

**W848**

Flat End Taper

3



Shank	L mm	ISO	Order No.	10,0
FG		• W806 314 173 534...	GW 848.314...	018

**W850**

Round End Taper

3



Shank	L mm	ISO	Order No.	10,0	10,0	10,0
FG		• W806 314 199 524...	W 850.314...	012	014	
		• W806 314 199 514...	FW 850.314...	012		
		• W806 314 199 534...	GW 850.314...	012	014	016

**W855**

Round End Taper

3




Shank	L mm	ISO	Order No.	7,0
FG		• W806 314 197 534...	GW 855.314...	025

**W856**

Round End Taper

3



Shank	L mm	ISO	Order No.	8,0	8,0	8,0	8,0
FG		• W806 314 198 534...	GW 856.314...	014	016	018	021

**W856L**

Long Round End Taper

3



Shank	L mm	ISO	Order No.	9,0
FG		• W806 314 199 534...	GW 856L.314...	018

**W858**

Needle

3



Shank	L mm	ISO	Order No.	8,0
FG		• W806 314 165 534...	GW 858.314...	014

**W860**

Flame

3



Shank	L mm	ISO	Order No.	5,0
FG		• W806 314 247 514...	FW 860.314...	012

**W862**

Flame

3



Shank	L mm	ISO	Order No.	8,0
FG		• W806 314 249 534...	GW 862.314...	014

**W863**


Flame

3




Shank	L mm	ISO	Order No.	10,0	10,0
FG		• W806 314 250 524...	W 863.314...	012	016
		• W806 314 250 514...	FW 863.314...	012	016
		• W806 314 250 534...	GW 863.314...	012	016

**W877K**  
Torpedo tapered  
3




Shank	L mm	ISO	Order No.	6,0	6,0
FG		• W806 314 297 524...	W 877K.314...	012	016
		• W806 314 297 514...	FW 877K.314...	012	016
		• W806 314 297 534...	GW 877K.314...	012	016

**W878**  
Cylindrical Torpedo  
3




Shank	L mm	ISO	Order No.	8,0	8,0	8,0
FG		• W806 314 289 524...	W 878.314...	010	012	014
		• W806 314 289 514...	FW 878.314...	010	012	014
		• W806 314 289 534...	GW 878.314...		012	014

**W878K**  
Torpedo tapered  
3




Shank	L mm	ISO	Order No.	8,0	8,0	8,0	8,0	8,0
FG		• W806 314 298 524...	W 878K.314...	014	016	018	020	022
		• W806 314 298 514...	FW 878K.314...	014	016	018	020	022
		• W806 314 298 534...	GW 878K.314...	014	016	018	020	022

**W879**  
Cylindrical Torpedo  
3




Shank	L mm	ISO	Order No.	10,0
FG		• W806 314 290 524...	W 879.314...	014
		• W806 314 290 514...	FW 879.314...	014
		• W806 314 290 534...	GW 879.314...	014

**W879K**  
Torpedo tapered  
3




Shank	L mm	ISO	Order No.	10,0	10,0	10,0	10,0
FG		• W806 314 299 524...	W 879K.314...	012	014	016	018
		• W806 314 299 514...	FW 879K.314...	012	014	016	018
		• W806 314 299 534...	GW 879K.314...	012	014	016	018

**W880**  
Round End Cylinder  
3




Shank	L mm	ISO	Order No.	6,0
FG		• W806 314 140 524...	W 880.314...	012
		• W806 314 140 514...	FW 880.314...	012
		• W806 314 140 534...	GW 880.314...	012

**W881**  
Round End Cylinder  
3



Shank	L mm	ISO	Order No.	8,0	8,0
FG		• W806 314 141 524...	W 881.314...	014	
		• W806 314 141 514...	FW 881.314...	014	016
		• W806 314 141 534...	GW 881.314...	014	016

**W899**  
Pointed Bull  
3



Shank	L mm	ISO	Order No.	7,0
FG		• W806 314 033 524...	W 899.314...	027
		• W806 314 033 514...	FW 899.314...	027
		• W806 314 033 534...	GW 899.314...	027

**W909**  
Round Wheel  
3



Shank	L mm	ISO	Order No.	1,6
FG		• W806 314 068 534...	GW 909.314...	040

