

GUTTA PERCHA OBTURATION SYSTEM



CV-FILL

Features precisely controlled motorized gutta percha extrusion, fast heating bendable thin tip, and digital control panel. Save valuable time with hygienic single-use gutta percha cartridge that provides reliable result every time. CV-FILL effectively and tightly compacts and seals all canals including lateral canals. It features LED light-guided condenser, instant heating, and digital temperature control. Both units can be fully charged within 90 minutes.

CV-FILL is the new advanced cordless obturation system that consists of CV-GUN obturation device and CV-PEN warm vertical compaction device. With the innovative electric motor,

GUTTA PERCHA OBTURATION SYSTEM

Feature:

Vertical Compaction Device obturation pen(CV-PEN)

- Quick heating tip reaches high level of temperature(200°C) within one second
- Excellent result with warm and vertical compaction for all canals including lateral canals. CV-PEN offers reliable results with tightly sealed root canals
- Automatically recognize tip size of heating condenser
- Adjustable Temperature Setting (140°C, 160°C, 180°C, 200°C) for safety
- Simple one-button operation

Obturation gun(CV-GUN)

- Controlled, precise motorized Gutta Percha extrusion
- Fast heating 200° , Gun needle can be 360 rotation , more easy to get the fill position
- Hygienic disposable cartridges to prevent secondary infections
- 3 different Temperature Settings (160°C, 180°C, 200°C) to control extrusion speed
- Easier, motorized extrusion and less hand fatigue
- No mess, easy clean up
- Li-ion chargeable battery:3.7v 2600mah
- Quick heating system: Temperature reaches 200°C within 23s
- The ergonomic design enables the obturation with a minimum input of force
- Operating time:continuous for 1.5h
- Gun needle:02,04,06


Package

- Order No : OKIT01 Regular Kit
 - CV-PEN Handpiece, CV-GUN Handpiece & Charger
 - Accessories: Injection Needle 3pcs, Heating Tips #F/FM/M Power Adapter & Cord, Scrubbing Brush
- *Gutta Percha Cartridge is not included in the Regular Kit.


Accessories

- Order No: Open01 CV-PEN Handpiece
- Order No: Ogun01 CV-GUN Handpiece
- Order No: OSP001 Injection Needle
- Order No: OSP002 Battery
- Order No: OSP003 Charger base
- Order No: OSP004 Power Adapter & Cord
- Order No: OSP005 Heating tips Kit (3 Pcs)
- Order No: OSP006 Heating Condenser #F
- Order No: OSP007 Heating Condenser #FM
- Order No: OSP008 Heating Condenser #M
- Order No: OSP009 Multi Tool
- Order No: OSP010 Scrubbing Brush

GUTTA PERCHA OBTURATION SYSTEM



GUTTA PERCHA OBTURATION SYSTEM



Radiographs



Perfect and tightly sealed root canls after warm obturation tecnique .