


Archwires NiTi Heat Activated

Pearson

Vivid NiTi Thermal (Pearson)

Specially designed heat activated Nickel Titanium wire which provides outstanding resiliency, true thermal performance, and consistent Af performance. Vivid Niti Thermal is soft at room temperature for easy ligation and is more resistant to permanent set. Centermark is indicated with a permanent etch mark. Premium Quality Made in USA



Round, Pkg of 10

Size	Upper	Lower
.012	U05-0280	U05-0282
.014	U05-0284	U05-0286
.016	U05-0288	U05-0290
.018	U05-0292	U05-0294
.020	U05-0296	U05-0298

Rectangular, Pkg of 10

Size	Upper	Lower
.014 x .025	U05-0300	U05-0302
.016 x .016	U05-0304	U05-0306
.016 x .022	U05-0308	U05-0310
.017 x .025	U05-0312	U05-0314
.018 x .018	U05-0320	U05-0322
.018 x .025	U05-0324	U05-0326
.019 x .025	U05-0328	U05-0330
.020 x .020	U05-0322	U05-0334
.021 x .025	U05-0336	U05-0338