


# Band & Bracket Removal

## Anterior Debonder/ Scaler

Royal Line

Our most popular bracket remover can also be used as an adhesive scraper. Heat-resistant Silicon tip is removable for corner-to-corner bracket removal with sharp wedge tips.

DE231



Tip D232

Dry Heat Sterilizable



## Diamond Sharpening File for Bracket Removers



DDS-1

Royal Line

## Posterior Band Remover Short

DE347SI



Tip D348TH

Dry Heat Sterilizable



Our most popular band remover has a medium-length sturdy jaw with a rounded tip. Great for use as a scraper, too.



Old-style Tip D348SH Dry Heat Sterilizable

Classic Line



D347I

## Posterior Band Remover Long

Long-style band remover with sharp-edged rounded removal tip.

NON-INSERTED



D347

## Posterior Band Remover Long

Non-inserted long band remover with pointed tip.

Royal Line



DE230N

## Angled Debonder Narrow

Angled version of our great wedge-edge bracket remover.

Classic Line



D230W

## Angled Debonder Wide

Wide jaws stay sharp longer. Great for occluso-gingival bracket removal.

Royal Line



DE160

## Anterior Debanding Plier

A superior instrument for scraping bonding material. Vinyl tip must be cold sterilized.

Tip D348LH

Dry Heat Sterilizable



Tip D161

