

## Anthos R7 Euro Swing

### Anthos R7 Chair

- > Memory foam option
- > Pneumatic headrest control option

### Anthos R7 Delivery system

- > Smart touch with digital display
- > Ceramic cuspidor bowl option
- > Doctor's side tray holder

### Venus Plus-LED Light

- > Head features 3-axis adjustment
- > Light intensity from 6,000 to 50,000 Lux
- > Color temperature 5,000K
- > Autoclavable handle

### Assistant's vac Pack

- > 1 3-Way Syringe
- > 1HVE , 1 SE



## Anthos R7 Standard Swing

### Anthos R7 Chair

- > Memory foam option
- > Pneumatic headrest option

### Anthos R7 Delivery system


- > Smart touch with digital display
- > Ceramic cuspidor bowl option
- > Doctor's side tray holder

### Venus Plus-LED Light

- > Head features 3-axis adjustment
- > On/off Switch
- > Light intensity from 6,000 to 50,000 Lux
- > Color temperature 5,000K
- > Autoclavable handle

### Rear Mount Vac Pack

- > 1 3-Way Syringe
- > 1HVE , 1SE





## Anthos R7 Mobile Cart System

### Anthos R7 Chair

- > Memory foam option
- > Pneumatic headrest control option

### R7 Cart

The Cart version ensures absolute working freedom. The mobility of the dentist's module is perfect for any workplace and allows optimal positioning of the operating instruments. Whether right-handed or lefthanded, the Cart solution offers an immediate response.

### Assistant's vac Pack

- > 1 3-Way Syringe
- > 1HVE , 1SE
- > Ceramic cuspidor option

### Venus Plus-LED Light

- > Head features 3-axis adjustment
- > On/off Switch
- > Light intensity from 6,000 to 50,000 Lux
- > Color temperature 5,000K
- > Autoclavable handle



## Anthos NEWBERG II Package

### Anthos R7 Chair

- > Memory foam option
- > Pneumatic headrest control option

### Newberg II Delivery system

- > The Newberg II combines design simplicity with style at a price that fits any budget. The Newberg II delivery system is a unique combination of flexibility and value. Available in post mount, swing mount or cart systems. If you are looking to change your delivery system and keep an existing chair then the NewbergII is the perfect choice.

### Forest 9045 Light

- > The Forest operatory light is a blend of style & function
- > The reflecting surface is made of neutral glass with a dichroic surface treatment
- > An integrated cooling fan keeps the light head cool to the touch

