



Galaxy Package

Nusimplicity Chair

- Thin tapered back
- Integrated touch pads located on both sides of upper backrest
- Easily visible illuminated "on" indicator
- Height from 15"-35.5" accommodates sit down and stand up dentistry.
- Optional airglide technology

Galaxy Delivery System

- Units available in fixed chair mounted or L/R version
- Piezo scaler, curing light and intraoral camera available integrated into delivery unit.

Everlight LED Light

- Natural daylight illumination provides clearer oral cavity visibility
- 9 color/intensity settings to meet all your operative needs, including color matching
- Consumes 70% less energy than halogens; lowers bills
- Available in track and fixed ceiling



Core Package

Core Chair

- Aluminum casting chair back, base, armrest, and headrest
- Double articulating headrest that adjusts easily to the patient

Core Delivery System

- Space for easy pressure monitoring
- Aluminum flex arm for stability and ease of positioning
- Optional intergrated dual - Space pads

Core LED Light

- Hands-Free adjustable light intensity to 30,000, 20,000 and 10,000 lux
- Programmable sensor for three brightness settings .

