

myQuickmat All-round

developed with **STYLE ITALIANO™**

new!

embrace
simplicity

Advantages

beautifully contoured
matrices for anatomical
restorations

comfortable matrix
positioning

stable matrix tightening

easy-to-tighten matrices

obstacle-free, accessible
operating field

ESPECIALLY INTENDED FOR RESTORATIONS IN CASE OF:

- Missing adjacent tooth
- Teeth with missing cusps
- Tooth isolation for core build-up
- MOD's

circumferential
matrix systems



polydentia
swiss manufacture 

polydentia.ch





new!

myQuickmat All-round

developed with **STYLE ITALIANO™**



embrace simplicity

Testimonial



"Modern conservative dentistry needs restorative solutions that simplify our procedures, making them more precise and efficient at the same time. myQuickmat All-round is a good example of this. It allows performing posterior restorations comfortably and straightforwardly, achieving anatomical results with natural contours and correct contact points."


Dr Jordi Manauta, Italy



Images courtesy of Dr Jordi Manauta




for more information



The Wedge

Interdental wedges that combine flexibility and rigidity for optimal matrix adaptation even in concave areas.



API is the Autoclavable Placement Instrument that makes Quickmat All-round stable and easy to place. Ergonomic and easy to use, API is the All-round Star!



myQuickmat All-round is the circumferential matrix system designed to perform anatomical posterior restorations of teeth with missing structures. The easy-to-use Quickmat All-round matrices enable restorations with naturally contoured anatomies and tight contact surfaces. Thanks to API (Autoclavable Placement Instrument), positioning Quickmat All-round matrices is stable and straightforward. API provides support and prevents the tensioning spool from tilting during matrix tightening. With myQuickmat All-round circumferential matrix system, you can comfortably treat challenging cases in the posterior region, working on an obstacle-free, accessible operating field.

kit



myQuickmat All-round introkit - REF 6850

- 2 pcs API (Autoclavable Placement Instrument for Quickmat All-round)
- 15 pcs Quickmat All-round circumferential posterior matrices, premolar - 0.03mm / 5.5mm
- 15 pcs Quickmat All-round circumferential posterior matrices, molar - 0.03mm / 6.3mm
- 45 pcs The Wedge interdental wedges - assorted

refills

REF 6851 API, Autoclavable Placement Instrument for Quickmat All-round 2 pcs	REF 6855 Quickmat All-round circumferential posterior matrices - premolar 0.03mm/5.5mm 35 pcs	REF 6856 Quickmat All-round circumferential posterior matrices - molar 0.03mm/6.3mm 35 pcs	
REF 7011 The Wedge Extrasmall (XS) 100 pcs	REF 7012 The Wedge Small (S) 100 pcs	REF 7013 The Wedge Medium (M) 100 pcs	REF 7014 The Wedge Large (L) 100 pcs