

PolydentiaPRO

new!

interproximal
tooth
protection
matrix

Advantages

excellent protection of
adjacent teeth during
preparation

contoured and 0.05mm thick

easy to insert and remove

less polishing time thanks to
scratch-free adjacent teeth



polydentia
swiss manufacture 

polydentia.ch



matrices

new!

PolydentiaPRO




interproximal tooth
protection matrix

Testimonial



“Certain work situations, in particular excavations or complex preparations, jeopardise the integrity of proximal surfaces. The protective matrix PolydentiaPRO is a ready-made, immediately available product suitable for these situations.”

Dr. Thomas Casanova,
Switzerland



Excellent protection of
the adjacent tooth against
inadvertent bur scratches




Clinical images by Dr. Giuseppe Chiodera



Easy insertion. Extremities can be used as a guide to pull the matrix.



Use a wedge to hold the matrix in place.



The adjacent tooth is now protected against bur scratches.

PolydentiaPRO is a pre-shaped, contoured protective matrix (0.05mm thick) designed to prevent damage to the adjacent tooth that commonly occurs during Class II preparations. The contoured shape and the two extremities are specially engineered to slide easily between the cervical margin and the contact point, facilitating matrix placement.



for more
information



PolydentiaPRO - REF 5720
Interproximal tooth protection matrices,
100 pcs

