


POLIDENT CAD-CAM MODEL DISCS



Suitable for milling dental models

Characteristics: This is a special, highly-compressed material with exceptional abrasion resistance. Recommended for ceramic layering. Nevertheless, it has excellent milling ability. Available in two densities.

Indications: Milling of stumps/ dies, working models and antagonists. Digital models are CAD designed from 3D data acquired with a wide range of digital impression-taking solutions.



TRANSPARENT CAD-CAM DISCS



Can be used for milling castable frameworks for metal casting and castable layers for press ceramic

Characteristics: Produced of 100% organic material that burns-out without any residual ash. Diameter: 95 - 105 mm. Thickness: 10-25 mm. Colour: transparent.

Indications: milling frameworks for casting, milling castable layers for press ceramic, verifying frameworks on plaster model before final ZrO₂ milling



POLIDENT d.o.o. Dental Products Industry
Vol ja Draga 42
5293 Vol ja Draga
Slovenia

Tel.: +386 5 330 48 40
Fax: +386 5 330 48 70
e-mail: polident@polident.si

www.polident.si

PoliDent



PoliDent
CAD CAM DISCS

High quality. Time saving. Esthetic.

POLIDENT PMMA CAD-CAM DISCS



Suitable for fabrication of temporary crowns and bridges and for verifying the bridge on the plaster model/ in mouth before the final ZrO₂ milling process

Characteristics: Vicker hardness: 26,60 HV. Flexural strength: 114 Mpa. E-Module: 2771 Mpa. Residual monomer < 1%. Diameter: 95 - 105 mm. Thickness: 10 - 25 mm.

Indications: temporary crowns and bridges, verification of the bridge on the plaster model/ mouth, splint

Colours: A1, A2, A3, B1, B2, B3, BI1, BI2, BI3, clear
enamel colours: E1, E2

Upon request we can supply also other shapes, diameters and/ or designs

POLIDENT PINK CAD-CAM DISCS



Can be used to mill digitally designed highly esthetic and functional dentures

Characteristics: Vicker hardness: 26,60 HV. Flexural strength: 114 Mpa. E-Module: 2771 Mpa. Residual monomer < 1%. Diameter: 95 - 105 mm. Thickness: 25 - 30 mm.

Indications: full dentures, partial dentures, temporary immediate load implant denture

Colours: pink

POLIDENT X-RAY CAD-CAM DISCS



X-ray visible disc used for milling diagnostic stents for implant planning

Characteristics: Contains X-ray visible powder. Homogeneous mix ensures best visibility. Prevents light scattering.

Indications: allows visibility of teeth setup in the X-ray CAT scan for implant planning and determination